

# ALL-ROUND TOXKIT INCUBATOR



**SUITABLE FOR ALL TOXKIT TESTS**

PERFORMANCE OF :

UP TO 4 (6 or 8) ALGALTOXKIT TESTS

UP TO 16 (24 or 32) TESTS (HATCHING & INCUBATION)

FOR ALL THE OTHER TOXKIT TESTS

## TECHNICAL DETAILS

<u>capacity</u> :	135 Ltr. (255 or 455 Ltr.)
<u>overall dimensions</u> :	600 deep x 600 wide x 850 high (mm) (600 x 610 x 1640 or 750 x 730 x 1640)
<u>inside dimensions</u> :	513 deep x 441 wide x 702 high (mm) (470 x 440 x 1452 or 600 x 560 x 1452)
<u>weight</u> :	42 kg (52 or 82 kg)
<u>construction</u> :	fully insulated cabinet with universal plug-in temperature
<u>control unit operation</u> :	robust, water protected foil front panel, 2 keys with tactile acknowledgement
<u>range of adjustment</u> :	2°C to 40°C; adjustable in 1°C-increments
<u>ambient temperature</u> :	10°C - 35°C
<u>display</u> :	3 digit for actual temperature, resolution 0,1°C 2 digit for temperature set-point
<u>refrigerating unit</u> :	compressor-system
<u>heating unit</u> :	heating element
<u>fan</u> :	radial, discharge capacity 120 m <sup>3</sup> /h
<u>power requirement</u> :	230V/50Hz
<u>light unit</u> :	4 (6 or 8) LED lamps (32.2 cm long) of 3.8 W that can be used individually
<u>power supply inside</u> :	3 (or 4) extra individual plug-sockets are available