

# STAND WITH TELESCOPIC ARM FOR A DIGITAL CAMERA TO TAKE PICTURES OF PHYTOTOXKIT TEST PLATES AND SPIRODELA TOXKIT MULTIWELLS

The robust stand is composed of 3 parts :

1. a horseshoe shaped base in cast iron, provided with a hinge to which the telescopic arm is fixed
2. a chrome-plated telescopic arm with a length range of 35-60 cm
3. an adjustable swivel shoe mount fixed at the end of the telescopic arm and provided with a bolt fit for screwing any type of digital camera

The hinge on the base of the stand and the swivel shoe mount allow to position the camera horizontally at any height above Phytotoxkit test plates and Spirodelela multiwells placed on a desk, for optimal shooting of pictures of the test containers.

