

PHYTOTOXKIT ACCESSORIES PACKAGE

This package is not included in the Phytotoxkit and has to be ordered separately.



The Phytotoxkit accessories package contains small items for easy test preparation, which can be used permanently.

The package contains the following items :

- a 100 ml graduated polypropylene beaker, for measurement of the 90 cm³ soil to be put into the test plates and for determination of the water saturation volume of the test soils
- a 50 ml graduated polypropylene cylinder to be used for determination of the water saturation volume of the test soils
- a 50 ml synthetic syringe, for hydration of the soil samples in the test plates
- a round (12 cm diameter) metal sieve of 2 mm mesh for sieving the test soils prior to use for the tests
- a spherical (5.5 cm diameter) wooden pestle to facilitate sieving of the test soil
- a small (3 cm diameter) flat pestle to compact sieved test soil in the beaker
- a metal spatula with a 4 cm blade in stainless steel and a handle, to spread the soil samples evenly in the lower compartment of the test plates and to flatten the soil surface.
- pair of plastic tweezers for easy transfer and positioning of the seeds on the paper filter in the test containers
- twenty rectangular plastic strips (2 cm height) to be placed in the bottom compartment of the test plates to facilitate the flattening of the soils
- a small pestle with flat bottom for flattening the soils in test plates provided with a rectangular strip